



Molift Hi-Trac Scale



NEW Molift Hi-Trac Scale

- **weight up to 300 kg / 47 st, accurate to 100 g**
- **calculates Body Mass Index (from 35 kg / 5st)**
- **battery and power saving - shuts down automatically if not used for 2 minutes**

molift[®]
designed for life

.. when accuracy counts.

- Can be used with Hi-Trac 2002/2004 (160 kg / 25st) and 2002+/2004+ (300 kg / 47st)
- Battery and power saving, the system will shut down if not used for 2 minutes
- Max weight 300 kg/47st, min weight 2 kg/4,4lbs
- The scale is to be fitted between the suspension and the lifting belt
- Easy to find the BMI - just press the BMI button until the display asks you to state the height in meters. Use the up/down arrows. Press the BMI button again to display.
- Easy to install, no tools required
- To be used with both slings and stretchers

Specification

Run by a replacable 9 V battery
type PP3

Displays the weighs in KG

Body Mass Index facility



Molift - designed for life

At times with increasing demand for efficiency in the healthcare sector and a growing focus on patient needs, it is easy to neglect the needs of the carer. The need to give better help and caring avoiding heavy lifting situations, avoiding putting your own back at risk, without being completely exhausted after work and still have some surplus energy when the workday is over.

This is the insight that drives us at Molift. The company which was started 20 years ago by a man who through his own experience saw the need for lifting and moving patients in an effective and comfortable way – for both patient and carer.

We have since then delivered products that compliment the natural movement of the body, that is both simple and intuitive to use. Products designed to give the patient a better life in addition to providing the carer with real support and optimum work enjoyment.

This is the true meaning of; Molift – designed for life

molift